
Pokyny	CZ
Bedienungsanleitung	DE
Brugsvejledning	DK
Instrucciones	ES
Käyttöohje	FI
Mode d'emploi	FR
Instructions	GB
Οδηγίες	GR
Használati utasítás	HU
Manuale d'Uso	IT
Gebruiksaanwijzing	NL
Bruksanvisning	NO
Instrukcja obsługi	PL
Instruções	PT
Bruksanvisning	SE
Uputstvo za upotrebu	YU

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Bellman Visit 868 Bridge, BE1520

Thank you for choosing products from Bellman & Symfon.

The Bellman Visit 868 System consists of a number of radio transmitters and receivers.

The transmitters detect different events in the surrounding area and transmit a radio signal to the receivers. The receivers pick up this signal and provide indications using light, sound and/or vibration.

The transmitter determines what type of light, sound or vibration should be displayed so that the reason for the indication is evident.

Read through the entire user manual first and then start to install the system.

Refer to the illustration of the Bellman Visit 868 System on the inside of the cover.

This is how the Bellman Visit 868 Bridge works

The Bellman Visit 868 Bridge conveys signals between the Bellman Visit System and the Bellman Visit 868 System. With a Bellman Visit 868 Bridge, you can therefore use products from the old and new system side by side.

Getting started

Unpacking, installing and testing the unit

1. Open the battery cover (7). Connect the power supply unit to the socket (6), fit one back-up battery, either a 6LR61 alkaline or 6F22 lithium type battery, and close the battery cover.
2. Press the Test Button (2). The Bellman Visit 868 Bridge lights up the LED (1) to indicate that the unit is transmitting a radio signal. The Bellman Visit 868 and Bellman Visit receivers will indicate the Door Signal.
3. Mount the Bellman Visit 868 Bridge on the wall either using the self-adhesive Velcro tape or suspended on the wallmount bracket (8) on the screw supplied. If the Velcro tape is used, the wall on which the transmitter is to be positioned can be cleaned using the enclosed wet wipe. The Bellman Visit 868 Bridge can also be placed on a level surface, e.g. a table, but you should avoid placing it on the floor.

Function

General

The BE1520 Bellman Visit 868 Bridge is a product for indoor use and operates as an interpreter, which conveys radio signals between the Bellman Visit and Bellman Visit 868 systems. With a Bellman Visit 868 Bridge, you can therefore use products from the old and new system side by side.

The Bellman Visit System will not of course have the same wide range of new indication options

as the new Bellman Visit 868 System, but the system will emit signals.

The Bellman Visit 868 Bridge also operates as a multi-function transmitter, either via a pushbutton or via an input for an Exterior External Trigger. The Pushbutton and the Exterior External Trigger can be detected separately and transmit different Signal Patterns to the receivers in the Bellman Visit 868 System depending on what has activated the Bellman Visit 868 Bridge. There is a range of options for connecting the Bellman Visit 868 Bridge to various applications.

Radio Key

On delivery all Bellman Visit and Bellman Visit 868 products are tuned to the same Radio Key. If you have a neighbour with a similar system, you can change to different Radio Keys so that you do not affect each other's systems. If you use the Radio Key Switch (5) to change the Radio Key on this transmitter, you must also change all other units in your Bellman Visit 868 System and your Bellman Visit System to the same Radio Key. Refer to the user manual for the relevant unit.

Indicators and Signals

System indicators

The LED (1) blinks, when the Bellman Visit 868 Bridge transmits or receives radio signals.

Power supply

When the Bellman Visit 868 Bridge is activated, the LED (1) normally blinks green. This means that the back-up battery is in good condition.

If the LED (1) blinks yellow, this means that the back-up battery is flat and must be changed. Only use a 6LR61 (alkaline) or 6F22 (lithium) type battery.

When the Bellman Visit 868 Bridge is connected to the power supply unit, the LED is constantly green, apart from when it is receiving or transmitting radio signals when it blinks.

Troubleshooting in brief

Problem	Solution
Nothing happens when the transmitter is activated with the test button.	<ul style="list-style-type: none"> • Check that the power supply unit is connected correctly. • If the power supply unit is not connected, change the back-up battery. Only use an alkaline 6LR61 or a lithium 6F22 type battery. <p>Make sure that the power supply unit is connected correctly.</p>
The LED (1) blinks yellow when the Bellman Visit 868 Bridge is activated!	<ul style="list-style-type: none"> • Change the back-up battery. Only use an alkaline 6LR61 or a lithium 6F22 type battery. <p>Make sure that the power supply unit is connected correctly.</p>
The LED (1) blinks green when the Bellman Visit 868 Bridge is activated but the receivers are not responding.	<ul style="list-style-type: none"> • Check the batteries in the receivers. • Check that the receivers are not placed too far away by moving them closer to the Bellman Visit 868 Bridge. • Check that the Bellman Visit 868 Bridge is set to the correct Radio Key. For further information see Function/Radio Key.

The receiver signals when no transmitter is activated.

- Change the Radio Key on all units in the system. For further information see Function/Radio Key.

1. LED. Combined transmission and back-up battery indicator. For further information see Indicators and Signals
2. Test button/pushbutton
3. Exterior External Trigger.
4. Signal Switch

5. Radio Key Switch
6. Connector for power supply unit
7. Battery cover
8. Wallmount bracket

For further information about the product in English, refer to the Appendix.

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Appendix - Further information

Connection

Connecting an Exterior External Trigger


The BE1520 Bellman Visit 868 Bridge has the following methods of activation:

- various accessories can be connected via the Exterior External Trigger socket (3), see below.
- with pushbutton

For more information about the indications provided by the receivers in the Bellman Visit 868 System, refer to Appendix/Further information/Settings below.

An Exterior External Trigger on the BE1520 Bellman Visit 868 Bridge can be connected to the vibrator output on several of Bellman & Symfon AB's products or to other equipment that produces a voltage or contact for activation. In addition, the following accessories can be used:

- magnetic switch BE9023 (accessory)
- door entering contact mat BE9024 (accessory)
- a standard pushbutton.

1. Entry monitor
2. Magnetic switch
3. Pushbutton
4. Door phone
5. 3.5 mm jack plug
6. Exterior External Trigger

For further information see Appendix/Further information/Technical information/Activation.

Settings

No adjustments are required for normal use. The relevant descriptions are provided below, if you wish to change a setting for some reason.

Radio Key

In order to use several Bellman Visit 868 and/or Bellman Visit Systems close to one another without interference, different Radio Keys can be set on the different systems. All Bellman Visit 868 System units are supplied from the factory tuned to the same Radio Key, channel 0. This means that all Radio Key Switches (5) on the transmitters are set to the OFF position. The Radio Key setting (referred to as a radio channel in the Bellman Visit System) is identical for both systems.

- To alter the Radio Key, move the Radio Key Switches (5) to the desired positions.

Please note: all Bellman Visit 868 and Bellman Visit products within a system must be tuned to the same Radio Key in order to operate as a group. However, all Bellman Visit receivers are activated when a fire alarm signal is detected regardless of the Radio Key.

Signal Pattern

A Signal Pattern is the name for the way in which a receiver in the Bellman Visit 868 System indicates activation. Changing the transmitters' Signal Switch changes the Signal Pattern which the receivers display when the transmitter is activated. Follow the instructions under Appendix/Further information/Settings/Setting the activation pattern to make your choice.

The following signal patterns for the Bellman Visit 868 System are available to the Bellman Visit 868 Bridge:

Type	LED-pattern	Sound	Vibration	Flash
Green 1	Green is constantly lit	1 x ding dong, lowfrequency tone	Separate	Yes
Green 2	Green blinks in sequences of two	2 x ding dong, lowfrequency tone	Separate	Yes
Green 3	Green blinks in sequences of three	1 x ding dong, highfrequency tone	Separate	Yes
Green 4	Green blinks constantly	2 x ding dong, highfrequency tone	Separate	Yes
Yellow 1	Yellow is constantly lit	1 x ring, lowfrequency tone	Short	Yes
Yellow 4	Yellow blinks constantly	2 x ring ring, highfrequency tone	Short	Yes
Orange 1	Orange is constantly lit	Baby	Rapid	Yes
Orange 2	Orange blinks in sequences of two	Baby	Rapid	Yes
Orange 3	Orange blinks in sequences of three	Baby	Rapid	Yes
Orange 4	Orange blinks constantly	Baby	Rapid	Yes
VMA	Red and Orange constantly blink alternately	VMA constant	Long	Yes
Fire alarm	Red blinks constantly	Fire alarm constant	Long	Yes

The following signal patterns for the Bellman Visit System are available to the Bellman Visit 868 Bridge:

Type	LED-pattern	Sound	Vibration	Flash
Green	Green is constantly lit	1 x ding dong, lowfrequency tone	Long	Yes
Yellow	Yellow is constantly lit	1 x ring, lowfrequency tone	Short	Yes
Orange	Orange is constantly lit	Baby	Long	Yes
VMA	Red and Orange constantly blink alternately	VMA constant	Constant	Yes
Fire alarm	Red blinks constantly	Fire alarm constant	Constant	Yes

Setting the Activation Pattern

The Bellman Visit 868 Bridge can be activated in two different ways.

On delivery the Bellman Visit 868 Bridge transmits the same signal regardless of the input causing the activation.

By changing the Signal Switches (4), different inputs on the same Bridge can transmit different Signal Patterns.

In this way the same Bridge can be used to forward signals from the Bellman Visit System to the Bellman Visit 868 System and vice versa at the same time as it can transmit a Door Signal when activated via the Exterior External Trigger (3) and a Baby Cry signal when activated via the Test Button (2).

The figure below shows the standard settings that are available:

Signal switch	Signal Pattern for the Bellman Visit 868 System when activated by:		Signal pattern for the Bellman Visit System when activated by:	
	Exterior External Trigger	Pushbutton	Exterior external trigger	Pushbutton
	Green 1	Green 1	Green	Green
	Green 2	Green 2	Green	Green

Signal switch	Signal Pattern for the Bellman Visit 868 System when activated by:		Signal pattern for the Bellman Visit System when activated by:	
	Exterier External Trigger	Pushbutton	Exterior external trigger	Pushbutton
	Green 3	Green 3	Green	Green
	Green 4	Green 4	Green	Green
	Orange 1	Green 4	Orange	Green
	Orange 2	Green 3	Orange	Green
	Orange3	Green 2	Orange	Green
	Orange 4	Green 1	Orange	Green
	Green 1	Orange 1	Green	Orange
	Green 2	Orange 2	Green	Orange

Signal switch	Signal Pattern for the Bellman Visit 868 System when activated by:		Signal pattern for the Bellman Visit System when activated by:	
	Exterier External Trigger	Pushbutton	Exterior external trigger	Pushbutton
	Green 3	Orange 3	Green	Orange
	Green 4	Orange 4	Green	Orange
	Green 1	Yellow 1	Green	Yellow
	Green 4	Yellow 4	Green	Yellow
	VMA	VMA	VMA	VMA
	Fire alarm	Fire alarm	Fire alarm	Fire alarm

Transferring signal patterns between systems

Since the Bellman Visit 868 System has more indication methods than the older Bellman Visit System, the systems will not give the same indications, for obvious reasons. The similarities will,

as far as possible, mean that the same LED colour will be used. The red LED is reserved for a fire alarm in the Bellman Visit 868 System, so that if the red LED is used by the Bellman Visit System, the indication given by the Bellman Visit 868 System will use a different colour. See the tables below.

The tables below show how the different Signal Patterns are transferred between the different systems.

From 868	To 433
Green 1	Green
Green 2	Green
Green 3	Green
Green 4	Green
Yellow 1	Yellow
Yellow 2	Yellow
Yellow 3	Yellow
Yellow 4	Yellow
Orange 1	Orange
Orange 2	Orange
Orange 3	Orange
Orange 4	Orange
VMA	VMA
Fire alarm	Fire alarm

From 433	To 868
Green	Green 1
Yellow	Yellow 1
Orange	Orange 1
Red	Green 4
Fire alarm	Fire alarm
VMA	VMA

Testing

It is easy to test the Bellman Visit 868 Bridge. If the Bellman Visit 868 Bridge does not work as described below, you can check further under Appendix/Further information/Troubleshooting/Troubleshooting guide.

How to test

A Bellman Visit 868 transmitter and receiver and a Bellman Visit transmitter and receiver are required to test the radio reception on the BE1520 Bellman Visit Bridge. All units must be tuned to the same Radio Key as the Bellman Visit 868 Bridge.

Testing the receiver in the Bellman Visit 868 Bridge:

- Make sure that the Bellman Visit 868 Bridge is connected correctly.
- Press the Bellman Visit 868 transmitter's test button
- The Bellman Visit 868 Bridge will indicate by blinking with the LED (1).
- Press the Bellman Visit transmitter's test button.
- The Bellman Visit 868 Bridge will indicate by blinking with the LED (1).

To test the Transmitter:

- Make sure that the Bellman Visit 868 Bridge is connected correctly.
- Activate via the Test Button (2) or via the Exterior External Trigger (3).
- The Bellman Visit 868 Bridge will now light up the LED (1) to show that it has been activated and that it is transmitting a signal to the receivers in the Bellman Visit 868 and Bellman Visit Systems.
- The receivers in both system will indicate an alarm according to the way that the Bellman Visit 868 Bridge has been set up using the Signal Switches (4).

Troubleshooting

You can carry out a number of checks yourself before sending a product for repair.

Troubleshooting guide

Problem	Solution
Nothing happens when the transmitter is activated with the test button.	<ul style="list-style-type: none"> • Check that the power supply unit is connected correctly. • If the power supply unit is not connected, change the back-up battery. Only use alkaline 6LR61 or lithium 6F22 type batteries. <p>Make sure that the power supply unit is connected correctly.</p> <ul style="list-style-type: none"> • Check that all the connections are correct.
The LED (1) blinks yellow when the Bellman Visit 868 Bridge is activated.	<ul style="list-style-type: none"> • Change the back-up battery. Only use an alkaline 6LR61 or lithium 6F22 type battery. <p>Make sure that the power supply unit is connected correctly.</p>

Problem	Solution
The LED (1) blinks green when the Bellman Visit 868 Bridge is activated but the receivers are not responding.	<ul style="list-style-type: none"> • Check the batteries in the receivers. • Check that the receivers are not placed too far away by moving them closer to the Bellman Visit 868 Bridge. • Check that all units in the Bellman Visit 868 and Bellman Visit Systems are set to the same Radio Key. For further information see Appendix/Further information/Settings/Radio Key.
One of the systems is working but the other one is not working.	<ul style="list-style-type: none"> • Check that all units in the Bellman Visit 868 and Bellman Visit Systems are set to the same radio key. For further information see Appendix/Further information/Settings/Radio key.
The receiver emits a strange alarm when activated by the Bellman Visit 868 Bridge.	<ul style="list-style-type: none"> • Set the Signal Switch to 0000.
The receiver signals when no transmitter is activated.	<ul style="list-style-type: none"> • Change the Radio Key on all units in the relevant Bellman Visit 868 System. There is probably another system nearby with the same Radio Key.

Technical information

Power supply

Mains power: 6 V DC / 800 mA with power supply unit BE9018 (Europe) or BE9083 (United Kingdom).

Back-up battery: 9 V 6LR61 alkaline

Power consumption: Active: 40 mA
Idle position: < 10 μ A

Radio function

Radio frequency: 868.3 MHz

Number of Radio Keys: 64 Radio Keys as standard. Special software can be used to increase these to 256 Radio Keys in increments of 64 per software purchase. Contact the nearest supplier for further information.

Coverage: For the Bellman Visit 868 system:
The normal coverage between a transmitter and receiver in the Bellman Visit 868 System is approximately 200 metres with a clear line of sight.

Coverage is reduced if walls and large objects screen off the signals. Any thick walls constructed of reinforced concrete will greatly affect coverage. The system may also be affected by radio transmitters such as TV transmitters, computers, mobile phones, etc. This means that a unit may



work perfectly in one part of the room but not at all in another.

For the Bellman Visit System:

The normal coverage between a transmitter and receiver in the Bellman Visit System is approximately 80 metres with a clear line of sight. Coverage is reduced if walls and large objects screen off the signal. Any thick walls constructed of reinforced concrete will greatly affect coverage.

Activation

Via test button

Exterior External Trigger: 3.5 mm stereo (mono provides a connection) jack plug (3)

Connection:

Between the inner and outer pins of the mono type 3.5 mm jack plug or between the middle/inner and outer pins of the stereo type 3.5 mm jack plug, see diagram.

DC: 2 to 30 V between the inner pin and middle pin on the stereo type 3.5 mm jack plug, see diagram.

AC: 3 to 24 V RMS 5 -150 Hz between the inner pin and middle pin on the stereo type 3,5 mm jack plug, see diagram.