### Inledning

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## Settings

Bellman Audio Mino has more settings than those you can access in user mode. By going into the Settings Menu you can adjust settings that you do not normally need to change or which are adjusted when first used and never need to be changed again.

### Settings Menu

How to get into the Settings Menu.


- Start the unit by pressing the On/Off button **(6)**.
- Hold down the On/Off button **(6)** and press down the Volume **—** button **(10)** within 1 second until the T LED **(4)** starts to blink.
- When the T LED **(4)** starts to blink, you can release both the buttons.
- Press and hold down the T button **(7)** within 2 seconds of the moment the T LED **(4)** starts to blink.
- Release the T button **(7)** immediately once the T LED **(4)** has gone out.

You are now in Balance Control mode.

### Overview of the Settings Menu

In the Settings Menu you can switch between the different settings by pressing the T button **(7)** or the MIC button **(5)**. Each press of the T button moves you clockwise on to the next setting in the list and the MIC button takes you anti-clockwise to the previous setting in the list.

When you have finished adjusting your settings, give a short press on the On/Off button **(6)** to return to normal user mode.

### Balance control

The balance between the right and left-hand channels can be set to compensate, for example, for a hearing loss in one ear.

Press the **<** button **(11)** to adjust the balance on the left and the **>** button **(12)** to adjust the balance on the right.

The standard setting is centered balance.

### Basic attenuation

Different headphones have different sensitivity levels and you can compensate for this with basic attenuation. What you actually set is the attenuation you require at the maximum power output you want to have (which is therefore 0 dB attenuation). This setting is useful when you want to adjust the maximum volume in user mode for a pair of headphones or in-ear phones with a higher sensitivity than the BE9122 Bellman Stereo Headphones or, for example, if you want to limit the volume for a child.

Use the **>** button **(12)** to reduce the basic attenuation (increase the volume) and the **<** button **(11)** to increase the basic attenuation (reduce the volume). The basic attenuation is adjusted in increments of 3 dB.

The standard setting is 0 dB attenuation, which corresponds to the maximum output signal.

### Stereo/Mono

The line-in signal can be rendered in either Stereo or in Mono. The intelligibility of speech can be increased in certain listening situations, such as when watching TV, if the signal is in Mono.

Press the **+** **(9)** button to select Stereo.

Press the **—** **(10)** button to select Mono.

The standard setting is Stereo.

### Bass reduction

Low-frequency interference such as fans and road noise can be tiring if you are exposed to it for a long period. By activating bass reduction, low-frequency interference and sounds are reduced.

Press the **+** **(9)** button to select Bass Reduction.

Press the **—** **(10)** button to deactivate Bass Reduction and get level frequency response.

The standard setting is deactivated Bass Reduction.

### Adjusting the Tone Control

The Tone Control Blocking function has been developed in case you want to make it easier for the user to use Mino. This setting can be used to choose whether or not the user should be able to adjust the Tone Control. If you choose to switch off the tone control buttons, you can preset the Tone Control that will be used in user mode.

### Activation and deactivation of Tone Control in User Mode

Press the **+** (9) button to activate the Tone Control buttons.

Press the **-** (10) button to deactivate them.

The standard setting is Tone Control activated.

### Presetting the Tone Control

Press the **<** button (11) and the **>** button (12) to preset the Tone Control to the level that will be used in user mode.

The standard setting is Tone Control centered (level frequency response)

*Please note that if the Tone Control is activated, the preset value will be changed as soon as the Tone Control is changed in user mode.*

### Blocking the T position

The T position can be blocked if you know that the T position will not be used. This setting is for those who will never use the T position and want to avoid selecting this position by mistake.

Press the **+** (9) button to activate the T position.

Press the **-** (10) button to deactivate the T position.

The standard setting is activated T position.

### Resetting to the factory setting

You can use this to reset all settings to the standard settings the product was delivered with from the factory.

Press and hold down the **+** button (9) until all the LEDs on the display (8) are lit. After resetting is completed, the display will blink for a short period.

### Version information

The product version information is presented here via a pattern of lit and blinking LEDs. This internal information is only intended as an aid for servicing and support.

## Testing

It is easy to test Mino yourself. If Mino is not working as described below, you can carry out further troubleshooting as instructed in the Troubleshooting section below.

### How to test

The testing of Mino has been divided into three parts. One part for internal microphones, one part for external sound and finally one for charging.

#### Internal microphones

- Switch on the unit and make sure that the Battery Status LED (3) lights up green.
- Connect the headphones to the Headphone Socket (13) and scrape carefully on the round microphone (1) with your finger and check that there is a scraping sound in the headphones.
- Press the **MIC** button (5) and make sure that the MIC LED (2) changes colour and lights up orange.
- Carefully scrape the small microphone hole on the top (1) and the round microphone (1) with your finger and check that there is a scraping sound in the headphones.
- Press the buttons (9, 10, 11, 12) and check that the Volume and Tone LED indicator changes.

#### External sound

- Press the **T** button (7) and make sure that the External Sound LED (4) lights up green.
- Hold the unit against an electrical watch with mechanical hands (with built-in battery) and check that it ticks.
- Connect an external sound source to the External Sound Source Input (16) and check that the sound from the sound source is heard in the headphones. The External Sound Source Input (16) can be tested, for example, by connecting a CD player or a radio to the input via a BE9126 audio cable. The unit automatically switches to an external sound source when a connector is connected to the External Sound Source Input (16).

#### Charging

- Ensure that Mino is switched off.
- Connect the charger to the mains socket.
- Connect the charging cable to both the charger and the unit's charging socket (14), see the fold-out section of the manual.
- Check that the unit's Battery Status LED lights up.

## Troubleshooting

Problem	Solution
The sound is louder in one ear than in the other.	<ul style="list-style-type: none"> <li>Adjust the Balance in the Settings Menu.</li> </ul>
The sound volume is low even though the volume is set to max.	<ul style="list-style-type: none"> <li>Adjust the Basic Attenuation in the Settings Menu.</li> </ul>
The sound is only in Mono even though it is connected via a stereo cable to a sound source which provides stereo sound.	<ul style="list-style-type: none"> <li>Change to Stereo in the Settings Menu.</li> </ul>
The bass reproduction is weak.	<ul style="list-style-type: none"> <li>Switch off Bass Reduction in the Settings Menu.</li> </ul>
The Tone Control cannot be adjusted.	<ul style="list-style-type: none"> <li>Use the Settings Menu to activate the Tone Control buttons or to adjust the sound.</li> </ul>
It is not possible to select T position.	<ul style="list-style-type: none"> <li>Use the Settings Menu to activate the T position.</li> </ul>
The unit is acting strangely.	<ul style="list-style-type: none"> <li>Reset the unit to the factory settings in the Settings Menu.</li> </ul>
Humming occurs when connecting the stereo cable while charging.	<ul style="list-style-type: none"> <li>Only use the supplied charger.</li> </ul>
Interference is heard when used in some environments	<ul style="list-style-type: none"> <li>Move the product away from the suspected source of radio frequency energy, for example a transmitter.</li> </ul>

## Technical information

Power supply	
Battery power:	Built-in rechargeable lithium-ion battery 3.7 V, 860 mAh
<p><b>CAUTION</b>  <i>Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type.</i></p>	
Charging device:	5V 1A DC adapter
Operating time:	BE9124 Earphones – up to 30 hours BE9159 Neck Loop – up to 10 hours
Charging time:	~ 4 hours
Input signals	
Microphone:	Built-in electric condenser microphones. (omnidirectional or directional)
Line input:	2.5 mm stereo jack plug, sensitivity 0 dBV (maximum input signal) Automatic selection between line input and external microphone
Telecoil:	Angled at 45° for optimal function horizontally and vertically
Output signals	
Headphone socket:	3.5 mm stereo jack plug
Output power headphone socket:	100 mW at 16 Ω load
Distortion:	
BE9122 Headphones	0.6% THD
BE9124 In-ear phones	0.08% THD
Dynamic range:	90 dB SNR (electrical)
Frequency range:	40-10,000 Hz (electrical)
Max. output level, saturation:	
BE9122 Headphones	124.3 dB SPL (@ 2.5 kHz)
BE9124 In-ear phones	132.0 dB SPL (@ 5 kHz)
BE9403 Stetoclip	145.0 dB SPL (@ 1kHz, 50 mW)

Amplification, full volume  
(1 kHz, 60 dB SPL input  
signal):

BE9122 Headphones	42.5 dB
BE9124 In-ear phones	55.0 dB
BE9403 Stetoclip	69.7 dB

Amplification, full volume  
(60 dB SPL input signal,  
HFA):

BE9122 Headphones	46.5 dB
BE9124 In-ear phones	51.5 dB
BE9403 Stetoclip	65.7 dB

### Properties

<b>Volume control:</b>	9 increments of 5 dB
<b>Tone control:</b>	+/- 10 dB in 9 increments, cut-off frequency 1.6 kHz
<b>On/Off button:</b>	Clear, countersunk pushbutton
<b>Microphone and T position:</b>	Separate pushbuttons
<b>Combined sound sources:</b>	For mixing external sound and sound from the internal microphone
<b>Balance control (right/left):</b>	0-20 dB programmable reduction in increments of 5 dB
<b>Setting for max. volume:</b>	0-21 dB programmable reduction in increments of 3 dB
<b>Bass Reduction setting:</b>	Lowering of frequencies under 400 Hz by 10 dB
<b>Dynamic compressor:</b>	<ul style="list-style-type: none"> <li>• 10 channels</li> <li>• Dynamic amplification 0-35 dB</li> <li>• Compression 3:1</li> <li>• Expansion 1:1.5</li> <li>• Attack time 20 ms &lt;3kHz, 25 ms ≥3kHz</li> <li>• Recovery time 110 ms &lt;3kHz, 200 ms ≥3kHz</li> </ul>

**Dynamic noise reduction:** 10 channel adaptive noise reduction

**Feedback handling:** Adaptive feedback reduction

**Digital signal processing:** 22,050 kHz sampling frequency  
16-bit resolution in stereo

**Additional information:**

- Stereo/Mono setting
- Resetting to the factory setting
- Casing designed to reduce contact noise

### Additional information

**Environmental requirements:** For indoor use only

**Operating temperature:** 0°–35° C, 32°–95° F

**Transport and storage temperature:** -10°–50° C

**Relative humidity:** 15%–90% (non condensing)

**External dimensions:** H × W × D 99 × 48 × 22 mm

**Weight:** 74 g

**Colour:** Grey with white front panel and silver-coloured buttons

### Regulatory requirements

FCC SDoC: FCC Part 15 Subpart B, CE, RoHS, WEEE, RCM

## Accessories

BE9122 Stereo headphones

BE9124 In-ear phones

BE9403 Stetoclip

BE9126 Stereo cable kit, 5 m (16 feet)

BE9127 External microphone, 5 m (16 feet)

BE9159 Neck loop


BE9129 MicroSet

BE9147 Travel case


## Measurements

Frequency range (acoustic) (BE9122 Headphones, BE9124 In-ear phones and BE9403 Stetoclip):


Preliminary measurements:

Output signal with BE9122,  
Stereo headphones

Output signal with BE9124,  
In-ear phones

Output signal with BE9403,  
Stetoclip

## FCC compliance statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Warning:** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

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
This Class B digital apparatus complies with Canadian ICES-003.

## Accessories

These original accessories are especially developed to maximize the performance of your Mino. Contact your nearest dealer or visit [bellman.com](http://bellman.com)

### Support centers

**Europe**  
Bellman & Symfon Group AB  
Södra Långebergsgatan 30  
436 32 Askim, Sweden  
Call +46 31 68 28 20  
E-mail [info@bellman.com](mailto:info@bellman.com)

**North America**  
Bellman & Symfon North America Inc.  
5509 Business Dr., Unit B  
Wilmington, NC 28405, USA  
Toll Free: (833) 235-5626

**Germany**  
Bellman & Symfon Deutschland GmbH  
Uelzener Str. 14  
21335 Lüneburg  
Call +49 800 7433338  
E-mail [info@bellman.de](mailto:info@bellman.de)

**Asia Pacific**  
Bellman & Symfon Asia Pacific  
1401-06 Main Tower, GD Int'l  
Building 339 Huanshi Dong Rd  
Guangzhou, China, 510095  
Call +86 20 8318 2950  
E-mail [info.cn@bellman.com](mailto:info.cn@bellman.com)

APP

## Warranty conditons

**Dear customer.** Thank you for your purchase of this Bellman & Symfon® product. Our vision is to supply outstanding products that are of great value for our customers. For this reason, we have put all our best efforts into designing and manufacturing this product. However, if you should have any questions about the product or if you require repair, please contact us.

**General warranty conditions.** Your Bellman & Symfon® product is covered by a limited warranty. Please read the following conditions carefully and keep this information in a safe place. Please also save your receipt or other proof of purchase as you will need this for warranty claims.

**Warranty period.** The warranty is depending on the type of product purchased and is valid from date of purchase.

Bellman & Symfon offers warranty as per below:

Bellman Products in general	Two years of warranty from date of purchase
Bellman Neckloop	One year of warranty from date of purchase
Rechargeable batteries	One year of warranty from date of purchase
Bellman Accessories	90 days of warranty from date of purchase

**Warranty conditions.** Bellman & Symfon warrants that the Product that you have purchased is free from manufacturing defects in materials and workmanship when dispatched from our warehouse. Your sales receipt is your proof of the date of purchase and the warranty period commences as of that date. During the warranty period, Bellman & Symfon will have the option for the defective product to be repaired or replaced. The replacement product will be of equal or greater function and in good working and physical order. Such a replacement will be supplied from refurbished stock. The repaired or replaced product will be under warranty for the remainder of the original product's warranty period.

**Warranty limitations.** This warranty automatically becomes void if the product has been physically damaged or rendered defective (a) as a result of accident, misuse, or abuse, or other circumstances beyond Bellman & Symfon's control; (b) by the use of parts or peripherals not authorized by Bellman & Symfon; (c) as a result of normal wear and tear; (d) by use within an improper operating environment; (e) by modification of the product; (f) by the serial number or product code sticker being removed or defaced; (g) as a result of service rendered by anyone other than a Bellman & Symfon Authorized Service Agent; (h) as a result of the product being operated other than in accordance with Bellman & Symfon Group AB's users manual. Neither does the warranty cover failure due to causes outside Bellman & Symfon's control such as failure in your phone carrier's line service (if applicable) or Acts of God. No warranty is offered for disposable batteries.

**Return policies / warranty claim.** Please see [www.bellman.com](http://www.bellman.com)