

Never miss a visitor, a call or
an important alarm again



RELIABLE ASSISTANCE WHEN YOU CAN'T DEPEND ON HEARING ALONE



DON'T LET YOUR HEARING LOSS HOLD YOU BACK

A hearing deficiency is like a fingerprint - unique to every person. Today's modern hearing aids are a great help to many, however they don't work for everyone. Even with a great hearing aid, there will be instances where they simply don't do the job. There are also situations where it is not practical to wear a hearing aid, for example when you are asleep.

Moreover, there are occasions when you need to relax and have some peace and quiet, without the fear of missing out on what is going on around you, for instance, if the telephone rings, the fire alarm goes off or if your baby demands love and support.

Bellman Visit will help you

notice important signals around you by transmitting them to where you are at that moment. They can also help you to use your other senses.

We have decided to call it the Bellman Visit Alerting system, a system that alerts you when you cannot rely on hearing alone.



AN EASY SOLUTION TO A DIFFICULT PROBLEM



10:18 COFFEE IN THE GARDEN
STAY REACHABLE ALL THE TIME

12:54 EXPECTING A CALL
ALWAYS KNOW WHEN THE PHONE RINGS

STAY IN TOUCH WITH YOUR FAMILY AND FRIENDS

Have you heard from friends or family that you are hard to reach over the phone or that you never seem to be at home when they come for a visit? Then it's time to have a look at what Bellman Visit has to offer.

Our Portable receiver and Pager receiver are great choices when you want to move around freely in your home, without running the risk of missing a call or a visitor. Just bring the receiver with you and you will be alerted with clear sounds, vibrations or flashing lights as soon as someone tries to reach or visit you. If you select a cord-

less receiver, you can even bring it outside.

All you need to do is to place the Door transmitter by your door bell and connect the Telephone transmitter to your telephone. That's all. No advanced programming or installation is required and you can start using the products immediately.



16:24 SURPRISE VISIT
WITH A DOOR TRANSMITTER YOU WILL NEVER MISS A VISITOR

▶ BE ALERTED INDOORS & OUTDOORS



Door transmitter



Telephone transmitter



Portable receiver



or Pager receiver

▶ BE ALERTED BY FLASHING LIGHT



Door transmitter



Telephone transmitter



Flash receiver

RELIABLE ASSISTANCE AROUND THE CLOCK



00:32

READY FOR A GOOD NIGHT'S SLEEP

NEVER OVERSLEEP AN IMPORTANT SIGNAL

BE RELAXED

WITHOUT BEING LEFT OUT

After a stressful day at work, it's nice to have some peace and quiet. If you are using a hearing aid and have to spend all day at work in a noisy environment, you will most likely remove it when you get home. With Bellman Visit, you can relax and enjoy the silence without being disconnected from the outside world.

Many people with hearing loss find it a bit worrying to go to bed because they fear missing a fire alarm or a baby cry. Our Alarm clock provides you with a safe and relaxing sleep. Aside from waking you up in the morning with your own choice of sound, vibration and flashing light, it also gives

you the comfort of knowing that you will never again miss an important signal during the night. A simple way of giving you and your loved ones a good night's rest.

All you need to do is to place a few of our Smoke alarms around the house and perhaps a Baby cry transmitter in the nursery.



04:24

FEEL SAFE AT NIGHT

KNOWING THAT YOU WILL WAKE UP WHEN YOU ARE NEEDED

> STAY SAFE FROM FIRE



Smoke Alarm transmitter




Alarm Clock Receiver



Bed shaker

> BE THERE FOR YOUR CHILD



Contact mat



Baby cry transmitter



Alarm clock receiver



Bed shaker

100% WIRELESS
NO INSTALLATION, EASY TO USE



WHY YOU SHOULD ALWAYS CHOOSE BELLMAN® VISIT

Bellman Visit is a totally wireless alerting system that gives you the freedom to move around in your house and garden as you please (approximate range of 30 meters). Thanks to the backup batteries, the system even works when there is a power failure.

Besides being totally wireless, Bellman Visit is a very easy and user friendly system. No installation is required. Just place the transmitter where you need it and bring a portable receiver with you when you move around. The Swedish quality and award winning design also give you a great

value for your money.

Bellman Visit is the simple wireless alerting system:

- ✓ Affordable
- ✓ Totally wireless
- ✓ No installation
- ✓ Easy to use
- ✓ Easy to carry
- ✓ Good quality
- ✓ Nice design



BELLMAN® VISIT WIRELESS ALERTING SYSTEM

TRANSMITTERS



TELEPHONE TRANSMITTER

Alerts you when the telephone rings. No installation required. Just connect it to the telephone socket.

BE1123



DOOR TRANSMITTER

Alerts you when someone rings your door bell. No installation required. Just place it close to your door bell.

BE1023



BABY CRY / VOICE TRANSMITTER

Alerts you when your baby or dependant needs you. No installation required. Just place it close to your baby or dependant.

BE1290




SMOKE ALARM TRANSMITTER

Alerts you if there is a fire. It is activated by smoke and sends a signal immediately to your receivers. (Appearance may differ from image.)

*BE1282 (USA) / BE1288 (CANADA) /
BE1280, BE1285 (AUSTRALIA, NEW ZEALAND)*


ACCESSORIES



MAGNETIC SWITCH

Detects when a door or window is opened. Connect it to the Telephone transmitter

BE9023



CONTACT MAT

Detects when somebody steps on it. Connect it to the Telephone transmitter.

BE9024

RECEIVERS



ALARM CLOCK RECEIVER

The Alarm Clock receiver alerts you with a flashing light, a sound and a vibration (if a bed shaker is connected) when it's time to wake up. It also indicates, with lighted symbols, when the phone rings or when your baby needs you.

BE1345



PAGER RECEIVER

The Pager receiver alerts you with clear vibrations and lighted symbols. It is small, light weight and can easily be carried in your pocket or on your belt. The charger (BE1260) is an accessory.

BE1230



PORTABLE RECEIVER

The Portable receiver alerts you with amplified sounds and a flashing light. It has an adjustable volume up to 90 dB and is very easy to carry.

BE1033



FLASH RECEIVER

The Flash receiver alerts you with a flashing light and coloured symbols. You can easily attach a Bed shaker to it.

BE1330

ACCESSORIES



BED SHAKER


The Bed shaker alerts you through easily recognized vibrations for different signals and alarms. It can be placed for example in the bed or the sofa. It can be connected to all the receivers mentioned above.

BE1270



Swedish hearing solutions for everyday life

For 20 years Bellman & Symfon of Sweden has worked hard to improve the quality of life of people with hearing deficiencies. For us it is a commitment not only to develop solutions based on the latest digital technology, but also to include unique and user friendly features that make a difference in everyday life situations. Practical and innovative solutions that help people to overcome some of the daily challenges they face due to their hearing loss.



LOGO Distributor

Distributor:
Adress, Address
Address
Phone Fax
E-mail www-adress



Bellman & Symfon®

Head Office

Bellman & Symfon Europe AB
S. Långebergsgatan 30
421 32 Västra Frölunda
Sweden
Phone: +46 31 68 28 20
Fax: +46 31 68 28 90
E-mail info@bellman.com
www.bellman.com

China Office

Bellman & Symfon Asia Ltd.
1108, Main Tower, GD Int'l Hotel
339 Huanshi Dong Rd
Guangzhou, China, 510098
Phone: +86 20 8318 2950
Fax: +86 20 8333 5528
E-mail info@bellman.com
www.bellman.com