

Dedicated to your hearing

Bellman & Symfon of Sweden is dedicated to improve the quality of life for people with hearing loss. Our people and partners are devoted to this mission and we work closely with hearing healthcare professionals and leading experts to protect your hearing and make better hearing possible for everyone.

Bellman & Symfon
Södra Långebergsgatan 30
436 32 Askim
Sweden
Phone +46 31 68 28 20
E-mail info@bellman.com


© 2019, Bellman & Symfon AB, all rights reserved. Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB.

bellman.com

BE1221_001ART001_EN

Vibio

Bluetooth® bed shaker



A great way to start your day

Whether you are ready to drop after the daily grind or just need a catnap to recharge, the Vibio Bluetooth® bed shaker makes sure you get up on time.

Set your alarm in the easy-to-use app and pop Vibio under your pillow or mattress. When it's time to start your day, it will wake you with silent but powerful vibrations.

Vibio is the perfect choice for heavy sleepers, couples and people with hearing loss.



Soft and cozy

Vibio features a soft and organic design with a quilted pattern that harmonizes seamlessly with your bed.

Easy to use app

Vibio connects to your mobile via Bluetooth® and comes with a free and user-friendly app that is available for both iOS and Android. In the Vibio app, you can create multiple alarms, choose vibration power and pick your favourite wakeup sound.

Multiple snooze options

Not ready to get out of bed yet? Snooze from the app or pull the fabric snooze strap under your pillow or mattress for a few more minutes in the sack.



Vibio highlights



Powerful vibrations
Wakes even heavy sleepers.



Bluetooth® 5
Reliable wireless communication.



Easy to use
User-friendly app for iOS and Android.



Notifications
Vibrates for incoming calls and messages.



Set and forget
Works even without the mobile.

